

Product Information

NerveRed C2 (FM4-64)

Catalog Number: N4073

Product Size: 5 mg, 5 × 1 mg

Application: Nerve terminal staining

Parameters

Appearance: Purple solid soluble in H₂O

Ex/Em (in MeOH): 543/- nm, (emission in MeOH is too weak to measure)

Ex/Em (in membranes): 510/750 nm

Molecular Formula: C₃₀H₄₅Br₂N₃

Molecular Weight: 607.51




Figure 1. Excitation and emission spectra of NerveRed C2 in liposomes.

Storage

Store at 4 °C and protect from light. When stored as directed, product is stable for at least 12 months.

Description

Nerve terminal probes are a series of cationic styryl fluorescent dyes used to track synaptic activity at neuromuscular junctions or synapses. These dyes usually have a lipophilic tail (two carbon chains) and a highly hydrophilic head with cations. The NerveRed probe is a dye with three double bonds.

Cationic styryl dyes function by staining synaptic vesicles.

When the dye is co-incubated with cells or tissues, the aqueous

part of the dye is not fluorescent, and the lipophilic tail of the dye is inserted into the cell membrane and shows strong fluorescence. After the nerve is stimulated, it will undergo endocytosis. At this time, the dye is wrapped in the vesicle, and the intensity of the fluorescent signal indicates the number of newly formed vesicles. Conversely, during exocytosis, the dye is released from the vesicle along with the neurotransmitter, resulting in a decrease in the fluorescent signal. Therefore, changes in fluorescence intensity reflect the situation of endocytosis / exocytosis or synaptic activity. The rate of fluorescence increase (binding rate) during endocytosis and the rate of fluorescence decrease (dissociation rate) during exocytosis vary depending on the type of dye. Generally, dyes with longer lipophilic tails and more double bonds have higher affinity for the membrane and therefore have higher binding rates and lower off-rates.

NerveRed C2 dye cannot be used for fixed cell staining. NM series dyes can be used for fixed cell staining.

Protocol

The following is a nerve terminal staining protocol for neuronal cells cultured on coverslips.

Nerve terminal dyes can also be used to label endocytic vesicles of non-neuronal cell types. The plasma membrane can be selectively labeled at 4°C, and the labeled endocytosis can occur within 10 min at room temperature or 37 °C. Tyrode





solution or other buffers can be used. The sodium ion channel blocker tetrodotoxin (TTX) can be optionally added to block the action potential and prevent the release of synaptic vesicles after staining. This protocol was optimized according to different experiments.

1. Dilute nerve terminal dyes to a final concentration of 4 μ M in 50 mM Tyrode solution. Place the coverslip containing the cells in this solution at room temperature for 1 min to completely submerge the cells.
2. Transfer coverslips to Tyrode + 0.5 μ M tetrodotoxin (TTX)

solution and incubate for 1 min at room temperature.

3. Wash the coverslips repeated at room temperature with Tyrode + 0.5 μ M TTX solution.
4. Observe with a fluorescence microscope.

Notes

1. There are quenching problems with fluorescent dyes. Please avoid light to slow down the fluorescence quenching.
2. For your safety and health, please wear lab coats and disposable gloves.

