

# Product Information

## Biotin cadaverine, trifluoroacetate salt

Catalog Number: B5034

Product Size: 10 mg

Application Scope: Carboxylic acid reactive dye


### Parameters

Appearance: White powder soluble in DMSO

Molecular Formula: C<sub>17</sub>H<sub>29</sub>N<sub>4</sub>O<sub>4</sub>S

Molecular Weight: 442.5

Molecular Structural:



### Storage

Store at 4°C and protect from light. When stored as directed, product is stable for at least 12 months.

### Description

Biotin cadaverine, trifluoroacetate salt can couple carboxylic acids, DNA and other biological molecules.

### Note

For your safety and health, please wear lab coats and disposable gloves.

