

Product Information

Biotin ethylenediamine iodoacetamide

Catalog Number: B5044

Product Size: 25 mg

Application Scope: Carboxyl reactive dye


Parameters

Appearance: White solid soluble in DMF

Molecular Formula: $C_{14}H_{23}IN_4O_3S$

Molecular Weight: 454.33

Molecular Structural:



Storage

Store at -20°C and protect from light. When stored as directed, product is stable for at least 12 months.

Description

Biotin ethylenediamine iodoacetamide is easy to couple with mercaptan.

Note

For your safety and health, please wear lab coats and disposable gloves.

