

Product Information

6 × DNA Loading Buffer

Catalog Number: D2041

Product Size: 5 × 1 mL

Storage: Store at 4°C

Description

6 × DNA Loading Buffer is used for loading DNA Markers and sample on agarose or polyacrylamide gel. In the electrophoretic process, two indicators of Bromphenol blue and patent blue VF show that the migration direction of these two indicators was the same as that of nucleic acid, which is a rough indicator of electrophoretic process.

Protocol

Add 1/6 volume of 6 × DNA Loading Buffer to your DNA sample with a final concentration of 1 ×. After mixing, it could be used directly in agarose gel or polyacrylamide gel electrophoresis.

